


Use large block display to compare Bur identification system when using the schematic below for proper Bur usage

● VT5 Set    ○ VT8 Set    ● VS8 Set

| Densifying Mode CCW (800-1500) RPMs / Cutting Mode CW (800-1500) RPMs |         |         |                                |               |               |               |       |       |                           |  |               |               |               |               |               |       |                           |  |
|---|---------|---------|--------------------------------|---------------|---------------|---------------|-------|-------|---------------------------|--|---------------|---------------|---------------|---------------|---------------|-------|---------------------------|--|
| BioHorizons®  |         |         | Tapered Internal, Tapered Plus |               |               |               |       |       |                           |  |               |               |               |               |               |       |                           |  |
|   |         |         | Soft Bone                      |               |               |               |       |       |                           | Hard Bone (Mandible)   |               |               |               |               |               |       |                           |  |
|   |         |         |                                |               |               |               |       |       |                           | In densifying mode make sure your osteotomy is 1.0 mm deeper than the actual implant final length. In extreme hard bone, utilize DAC (Densify After Cut) Protocol. <a href="#">Click here</a> to view PDF. |               |               |               |               |               |       |                           |  |
| Geometry  | Major Ø | Minor Ø | Pilot                          | Bur 1         | Bur 2         | Bur 3         | Bur 4 | Bur 5 | Densah® Bur Block Display | Pilot  | Bur 1         | Bur 2         | Bur 3         | Bur 4         | Bur 5         | Bur 6 | Densah® Bur Block Display |  |
| Taper   | 3.0     | 2.0     | Pilot                          | VT1828* (2.3) | —             | —             | —     | —     |                           | Pilot  | VT1828* (2.3) | —             | —             | —             | —             | —     | —                         |  |
| Taper   | 3.0     | 2.8     | Pilot                          | VT1828 (2.3)  | VS2228* (2.5) | —             | —     | —     |                           | Pilot  | VT1828 (2.3)  | VS2228* (2.5) | —             | —             | —             | —     | —                         |  |
| Taper   | 3.5     | 3.1     | Pilot                          | VT1525 (2.0)  | VS2228* (2.5) | —             | —     | —     |                           | Pilot  | VT1525 (2.0)  | VT1828 (2.3)  | VT2535* (3.0) | —             | —             | —     | —                         |  |
| Taper   | 4.5     | 3.9     | Pilot                          | VT1525 (2.0)  | VT2535 (3.0)  | VS3238* (3.5) | —     | —     |                           | Pilot  | VT1525 (2.0)  | VT1828 (2.3)  | VT2535 (3.0)  | VT2838* (3.3) | VT3545* (4.0) | —     | —                         |  |

\*Denotes implant placement.

NOTE: Surgeon preference overrides this suggestive protocol




Use large block display to compare Bur identification system when using the schematic below for proper Bur usage

● VT5 Set    ○ VT8 Set    ● VS8 Set

| Densifying Mode CCW (800-1500) RPMs / Cutting Mode CW (800-1500) RPMs |         |         |   |              |               |               |               |       |                           |  |              |              |              |               |               |               |                           |
|---|---------|---------|---|--------------|---------------|---------------|---------------|-------|---------------------------|--|--------------|--------------|--------------|---------------|---------------|---------------|---------------------------|
| BioHorizons®  |         |         | Tapered Internal, Internal Tissue Level |              |               |               |               |       |                           |  |              |              |              |               |               |               |                           |
|   |         |         | Soft Bone                               |              |               |               |               |       |                           | Hard Bone (Mandible)   |              |              |              |               |               |               |                           |
|   |         |         |   |              |               |               |               |       |                           | In densifying mode make sure your osteotomy is 1.0 mm deeper than the actual implant final length. In extreme hard bone, utilize DAC (Densify After Cut) Protocol. <a href="#">Click here</a> to view PDF. |              |              |              |               |               |               |                           |
| Geometry  | Major Ø | Minor Ø | Pilot                                   | Bur 1        | Bur 2         | Bur 3         | Bur 4         | Bur 5 | Densah® Bur Block Display | Pilot  | Bur 1        | Bur 2        | Bur 3        | Bur 4         | Bur 5         | Bur 6         | Densah® Bur Block Display |
| Taper   | 3.8     | 2.8     | Pilot                                   | VT1525 (2.0) | VT2535* (3.0) | —             | —             | —     |                           | Pilot  | VT1525 (2.0) | VT1828 (2.3) | VT2535 (3.0) | VT2838* (3.3) | —             | —             |                           |
| Taper   | 4.6     | 3.9     | Pilot                                   | VT1525 (2.0) | VT2535 (3.0)  | VT3545* (4.0) | —             | —     |                           | Pilot  | VT1525 (2.0) | VT1828 (2.3) | VT2535 (3.0) | VT2838 (3.3)  | VT3545* (4.0) | —             |                           |
| Taper   | 5.8     | 3.9     | Pilot                                   | VT1525 (2.0) | VT2535 (3.0)  | VT3545 (4.0)  | VT3848* (4.3) | —     |                           | Pilot  | VT1525 (2.0) | VT2535 (3.0) | VT2838 (3.3) | VT3545 (4.0)  | VT3848 (4.3)  | VT4555* (5.0) |                           |

\*Denotes implant placement.

NOTE: Surgeon preference overrules this suggestive protocol



Use large block display to compare Bur identification system when using the schematic below for proper Bur usage

● VT5 Set ○ VT8 Set ● VS8 Set

| Densifying Mode CCW (800-1500) RPMs / Cutting Mode CW (800-1500) RPMs |         |         |           |              |               |              |               |               |                           |  |              |              |              |               |               |              |               |               |                           |
|---|---------|---------|-----------|--------------|---------------|--------------|---------------|---------------|---------------------------|--|--------------|--------------|--------------|---------------|---------------|--------------|---------------|---------------|---------------------------|
| BioHorizons®  |         |         | Internal  |              |               |              |               |               |                           |  |              |              |              |               |               |              |               |               |                           |
|   |         |         | Soft Bone |              |               |              |               |               |                           | Hard Bone (Mandible)   |              |              |              |               |               |              |               |               |                           |
|   |         |         |           |              |               |              |               |               |                           | In densifying mode make sure your osteotomy is 1.0 mm deeper than the actual implant final length. In extreme hard bone, utilize DAC (Densify After Cut) Protocol. <a href="#">Click here</a> to view PDF. |              |              |              |               |               |              |               |               |                           |
| Geometry  | Major Ø | Minor Ø | Pilot     | Bur 1        | Bur 2         | Bur 3        | Bur 4         | Bur 5         | Densah® Bur Block Display | Pilot  | Bur 1        | Bur 2        | Bur 3        | Bur 4         | Bur 5         | Bur 6        | Bur 7         | Bur 8         | Densah® Bur Block Display |
| Straight  | 3.5     | 3.5     | Pilot     | VT1525 (2.0) | VT2535* (3.0) | —            | —             | —             |                           | Pilot  | VT1525 (2.0) | VT1828 (2.3) | VT2535 (3.0) | VT2838* (3.3) | —             | —            | —             | —             |                           |
| Straight  | 4.5     | 4.0     | Pilot     | VT1525 (2.0) | VT2535 (3.0)  | VT2838 (3.3) | VS3238* (3.5) | —             |                           | Pilot  | VT1525 (2.0) | VT1828 (2.3) | VT2535 (3.0) | VT2838 (3.3)  | VT3545* (4.0) | —            | —             | —             |                           |
| Straight  | 5.7     | 5.0     | Pilot     | VT1525 (2.0) | VT2535 (3.0)  | VT3545 (4.0) | VT3848 (4.3)  | VS4248* (4.5) |                           | Pilot  | VT1525 (2.0) | VT1828 (2.3) | VT2535 (3.0) | VT2838 (3.3)  | VT3545 (4.0)  | VT3848 (4.3) | VT4555* (5.0) | —             |                           |
| Straight  | 6.0     | 6.6     | Pilot     | VT1828 (2.3) | VT2838 (3.3)  | VT3848 (4.3) | VT4858* (5.3) | —             |                           | Pilot  | VT1828 (2.3) | VT2535 (3.0) | VT2838 (3.3) | VT3545 (4.0)  | VT3848 (4.3)  | VT4555 (5.0) | VT4858 (5.3)  | VS5258* (5.5) |                           |

\*Denotes implant placement.

NOTE: Surgeon preference overrides this suggestive protocol